

Laboratory work Station

SEMCO offers the standard in workbench design and construction. Component for component, SEMCO furniture work station has been designed to offer you the ultimate in durability, flexibility, and performance. Whether it's critical weld locations or any of the many unique SEMCO design features, you'll be impressed by SEMCO's commitment to detail. When you demand quality design, specify SEMCO.



HOW TO USE LABORATORY WORK STATION:

- Always keep work station clean and dry.
- In case of spillage of chemicals, Wipe off the work surface immediately, Delay may damage the affected work surface.
- Turn off all heating apparatus, gas valves, And water faucets when not in use.
- Solvents should be in flame proof cabinets only.
- Ensure that drawers are not over loaded.
- Do not use laboratory work bench as a stool.
- Do not use work bench as dining table.
- Always wear appropriate eye protection.
- Long-sleeved laboratory coat and shoes that cover the whole foot.
- Deep scratch on powder coated component invites rust on exposed portion, immediate touch-up will ensure long life of the component.



CRCA / GI



STAINLESS STEEL

NOTE: Technical specifications & Design are subject to change without prior notice, due to continuous up gradation.

Pedestal Assembly



NOTE: Technical specifications & Design are subject to change without prior notice, due to continuous up gradation.

Installations



NOTE : Technical specifications & Design are subject to change without prior notice, due to continuous up gradation.